

# SATDrape

Universal Sterile CT Scanner Drape that reduces cross-contamination

The worlds first Sterile CT scanner drape that can be used in Interventional and routine CT patients



CE U.S. FDA Listed

## Key Benefits



### Easy to USE

30 seconds to Install - can be USED on all CT scanners - Universal



### Absorbent and Sterile

Can hold upto 2 liters of fluid for 12 hours. Can be used for Routine and Interventional CT



### Reduce Cleaning Times


Increase workflow by reducing cleaning times after Trauma, ICU and COVID patients

## Technical Specifications

SKU Number	72233
Product Name	Sterile Disposable CT Drape
Product Category	Surgical Drape
Feature	Sterile Single-USE and Routine Multi-USE
Absorption and Saturation	2 liters on the surface without leakage at maximum   Volume and saturation for 12 hours 1 liter in the collection bag without leakage at maximum saturation for 12 hours
Material	Non-woven fabric exterior and underlay (surface touching CT scanner surface) is non-permeable
Dimension (mm)	Front Piece 135 x 160 cm   Connected back piece 42 x 45 cm
Weight (g)	Box dimension: 60 x 48 x 32 cm   Box Weight: 13.5 kg/carton (Net Weight 12.5kg/carton)
Sterilisation	Ethylene Oxide (EO)
Shelf Life	3 or 5 Years
Packing	1 drape per pouch   50 pcs per carton

## Application & Placement

Can be used on all CT scanner manufacturers - Both Human and Animal CT Scanners




Place the SATDrape before each interventional CT procedure and one per day (in routine CT) until it is soiled for maximum workflow efficiency

## Product Structure

### Outer Layer

Non-Woven Fabric



### Inner Layer

Non-permeable membrane

### Overflow Bag

Can Hold extra 1 Ltr



Discover SATDrape

✉ sales@satmed-health.com

🌐 <https://satmed-health.com>

Note: The information in this technical sheet is subject to change without prior notice as part of our continuous product improvement policy.